

Masks

Anti-Viral Community Snood

This tubular snood is designed to stretch over the head and to cover the mouth and nose. Through the provision of a single seam to the back and raw edges, whilst using a fabric that will not fray, our snood is designed to provide enhanced wearer comfort. The snood fabric is finished with an antiviral treatment to protect the treated article from contamination against viruses and bacteria -ViralOff®. The Silver ions in the finish, strongly inhibit bacterial growth and viral activity, by interfering with the DNA. This kills both bacteria and viruses by 99% on a treated material over the space of two hours (ISO18184:2019 (H3N2)). It is not intended to cure or prevent diseases. This snood can be used as a preventative barrier. For best performance and sustainability, wash less and only when needed. Washing is not necessary as it is 'self-sterilizing'. If the product is washed we recommend hand wash or spot clean. However the product can be washed at 30 degrees. Please note this snood is neither a medical device or an item of approved PPE and must be worn and used in conjunction with hygiene guidelines from the health authorities.



Key features

N/A

Garment details

Style code	FMK5
Fabric composition	80% Polyester / 20% Elastane
Fabric weight	170gsm
Fabric care	Hand wash or spot clean. Can be washed at 30°

Performance standards

N/A

Colours

Black

Size range

One Size

Fits